



AIR LIQUIDE IN THE HOSPITAL

EZ-OX™ Portable Oxygen System



Air Liquide Healthcare is the leader for innovative medical gas products, being the first to provide a portable oxygen solution to hospitals.

YOU NEED PORTABLE OXYGEN - FAST

WHEN EVERY SECOND IS PRECIOUS, YOU WANT:

- A quick and easy oxygen set-up for your patients
- A reliable oxygen supply to get them to the Emergency Room, or to transport them across the hospital
- To keep your patients stable in a location where there is no headwall unit

You're frustrated with the delays and difficulties in finding and assembling regulators and oxygen cylinders, while your patient needs to be in surgery or back in recovery.

AIR LIQUIDE HEALTHCARE HAS YOUR SOLUTION

EZ-OX™ is the innovative solution: a pressure regulator, flow meter and medical E-oxygen cylinder all integrated into a single lightweight unit.

With no assembly needed and an intuitive, easy-to-use design, EZ-OX™ is the reliable oxygen bridge which nearly anyone can operate with one short training session.

Let EZ-OX™ be your new portable oxygen solution, so that you can focus on the job in hand—getting patients healthy.

WORKING HARD SO YOU CAN BREATHE EASY™

EZ-OX™ Portable Oxygen System



EZ-OX™

ESTIMATED USAGE PERIOD

	2 lpm	6 lpm	10 lpm	15 lpm
FULL	5 hr	1 hr 45 min	1 hr	40 min
3/4 FULL	3 hr 45 min	1 hr 15 min	45 min	30 min
1/2 FULL	2 hr 30 min	45 min	30 min	20 min
1/4 FULL	1 hr 15 min	20 min	15 min	10 min

INTERESTED IN TRUE PORTABLE OXYGEN PERFORMANCE?

We offer **EZ-Trials™** where you can try out **EZ-OX™** units at your facility. Contact your local Air Liquide Healthcare account manager for details or call Air Liquide Healthcare direct at **1-888-889-2793**.

EZ-OX™ - CLINICAL BENEFITS

- Integrated regulator and flow meter saves time with a quick, easy setup and instant readiness.
- Easy-to-use, lightweight design reduces operator error.
- Unique, balanced regulator design reduces incidences of the cylinder falling over.
- Protective cap over the regulator and flow meter increases unit reliability.
- Certified for use with 1.5 Tesla MRIs.

EZ-OX™ - COST EFFECTIVENESS

- Eliminates costly regulator purchases and repairs.
- Reduces training costs and patient setup time.
- Simplicity enables operation by all hospital staff—not just clinical personnel.
- Reduces accidents through simple design and user friendliness.

EZ-OX™ - KEY SPECIFICATIONS

- Integrated regulator
- Indexed flow meter with continuous flow using a graduated knob with a range from 2 to 15 liters per minute (lpm)
- Easy to read content gauge
- Carrying handle
- Protective cap
- Hose barb connection
- Usage chart label with safety instructions.
- No maintenance or repair required (Air Liquide provides full service)
- Fits into all existing racks and carts.

THE AIR LIQUIDE DIFFERENCE

- Medical gases and services to the highest standards
- Customer service second to none
- A total commitment to the patient and those who care for them

www.US.AirLiquideHealthcare.com



NOTE: This brochure is intended for general information purposes only and is not intended as a representation or warranty of any kind, or as a statement of any terms or conditions of sale. The information herein is believed to be correct, but is not warranted for correctness or completeness, or for applicability to any particular customer or situation. EZ-OX and EZ-Trials are trademarks of Air Liquide Healthcare America Corporation; © 2010 AIR LIQUIDE HEALTHCARE AMERICA CORPORATION